

Title: Electrographic Flow Consistency Correlates with High-Density Bipolar Voltage Maps in Sinus Rhythm and Atrial Fibrillation

Authors: David Haines, MD, Steven Castellano, PhD, Melissa Kong, MD, Philipp Sommer, MD

ABSTRACT:

Background

In patients with atrial fibrillation (AF), left atrial scar revealed by low voltage areas on 3D electroanatomic mapping have been correlated with scar on late gadolinium enhanced cardiac magnetic resonance imaging, which is recognized as a surrogate for fibrosis. Electrographic Flow (EGF) mapping reconstructs global atrial wavefront propagations in near real-time. Electrographic flow consistency (EGFC) quantifies concordance of observed wavefront patterns and may provide insight into the nature of atrial substrate where diseased substrate shows more chaotic wavefront propagation and lower FC.

Objectives

Correlate EGFC with regional bipolar voltage values in patients during AF and sinus rhythm (SR).

Methods

FLOW EVAL-AF is a single-center, prospective study of 10 AF patients who underwent concomitant high-density bipolar voltage mapping with HD Grid and EGF mapping with 64-pole basket catheter. Biatrial voltage maps and 1-minute EGF maps from 2-3 standard basket positions per atria were made prior to any induction or ablation. Atrial surface area with voltage >0.5 mV was computed from 2D voltage map projections. EGF fields computed per 2-sec segment of recording. Overall EGFC per atrium computed as average of the modulus of individual EGF vectors, where vector length represents consistency of flow patterns in arbitrary units. We compared % atrial surface area with normal voltage v. mean EGFC across EGF recordings per atria in a given rhythm.

Results

Data from 8 patients (2 AF; 6 SR) was analyzed; 2 excluded for non-comparable voltage map acquisition parameters. Mean age 70 ± 11 years; 1 female; mean LA size 4.5 ± 0.35 cm. Surface area % with voltage >0.5 mV was 0.24 ± 0.08 in AF v. 0.36 ± 0.05 in SR. Mean EGFC in AF 0.74 ± 0.14 and in SR 1.0 ± 0.11 . Intraprocedural rhythm change occurred in 5 (3 induced into AF; 2 converted to SR post-PVI), and all had higher voltage values overall and higher mean EGFC in SR v. AF. Pooling data from all pre-procedure voltage maps and EGF maps revealed overall correlation between EGFC and bipolar voltage from both atria and in both rhythms ($r=0.67$, $p=4.8 \times 10^{-3}$).

Conclusions

Both EGFC and bipolar voltage values were higher in SR v. AF. EGFC correlated with voltage across both atria and in both SR and AF. Global EGFC can be obtained with 1-minute recordings in 2-3 basket positions per atria and may provide another measure of the degree of atrial myopathy similar to the information gained from high-density bipolar voltage mapping.

Figure:

A) Top left: High-density bipolar voltage map of the left atrium (LA) generated using an HD Grid catheter in a patient who presented in persistent AF and at top right, the corresponding baseline LA EGF map recorded using a 64-electrode basket catheter (position shown in 3D electroanatomic map). Bottom left: New high-density bipolar voltage map of the LA in the same patient but now post-PVI with conversion to SR and at bottom right, the corresponding LA EGF map recorded with the 64-electrode basket catheter in essentially the same position within in the LA (shown in 3D electroanatomic map) as the map recorded in AF. In the EGF maps, higher flow consistency is indicated by longer flow vectors shown in magenta and correlated with normal voltage values (purple areas on voltage maps).

B) Plot of baseline pre-ablation high voltage fraction v. mean electrographic flow consistency. The presenting rhythm of the datum is indicated by its color, and an overall best-fit linear trendline across all data is indicated in orange ($r=0.67$, $p=4.8 \times 10^{-3}$). There are 8 data points from LA and 8 from RA.

